

Axis-Shield Liquid Stable 2-Part Homocysteine Reagent

ASSAY PROTOCOL – BECKMAN DXC 600-800

User defined setup

Parameters				
Analyte	[Homocysteine]			
Reaction Type	[Kinetics 1]			
Units	[µmol/l]			
Precision	[X.XX]*			
React. Direction	[negative]			
Curve Type	[linear]			
Primary Wavelength	[340] nm	Secondary Wavelength	[380] nm	
Factor	[1.000]			
No. of Cals.	[2]			
Set Point:				
1	[0.000]	4	N/A	
2	[28.000]	5	N/A	
3	N/A	6	N/A	
Cal. Limit	[336]			
Parameters				
	Primary	Secondary	Tertiary	
	A	NONE	C	
Quantity	[200]		[20]	
Time of Injection		[-180]	[270]	
Sample Quantity	[20]			
Reading Windows				
Reaction	White	Initial	Reaction 1	Reaction 2
Start Reading	[-50]	[273]	[357]	
End Reading	[-10]	[288]	[632]	
Error Detection Limit				
	White	Reaction 1	Reaction 2	
ABS Lower Limit	[-1.500]	[-1.500]	[-1.500]	
ABS Upper Limit	[2.200]	[2.200]	[2.200]	
Lower Limit Rate	[2.200]	[2.200]	[2.200]	
Upper Limit Rate	[-1.500]	[-1.500]	[-1.500]	
Mean Deviation	[2.200]	[2.200]	[2.200]	
Substrate Consumption				
Initial Kinetic	[-99.999]			
Delta ABS	[2.200]			
Multipoint Span				
1-2	[-0.001]			
Rang res. Util.		Interval ORDAC		
Lower Limit	[0.000]	[0.000]		
Upper Limit	[80.000]	[99999.999]		

* - User Defined

Axis-Shield Diagnostics Ltd.

Luna Place
The Technology Park, Dundee
DD2 1XA, Scotland, UK
Tel: +44 (0) 1382 422000
Fax: +44 (0) 1382 422088
Email: homocysteine@axis-shield.com



www.homocysteine.org.uk
www.axis-shield.com