

Axis-Shield Liquid Stable 2-Part Homocysteine Reagent

ASSAY PROTOCOL – Roche cobas 6000 (c501) / cobas 8000 (c502)

Test: HCY*	Type: Ser/PI
ANALYZE	
Assay time/Point	[2 Point End]/[10][40][70][0][0]
Wave (2nd/Primary)	[376]/[340]
S. Vol (Normal)	[9.0]
Reagent (R1) T1	[135] [0] [000000]
Reagent (R2) T2	[0] [0] [000000]
Reagent (R3) T3	[15] [0] [000000]
Linearity Limit	[30][10][60][60]
Prozone Limit	[0] [Lower]
Abs. Limit	[0][0][0][0][0][0][Lower][0][0]
Cell Detergent	[Detergent 1]
CALIB	
Calibration Type	[Linear]
Point	[2]
Span	[2]
Weight	[0]
Auto Calibration	
2 Point	[7]
SD Limit	[0.1]
Duplicate Limit	[10][32000]
Sensitivity Limit	[-99999] [99999]
S1 Abs. Limit	[-32000] [32000]
RANGE	
Application Code*	[] Unit [µmol/L]
OTHERS	
<Standard>	(1) (2)
Calibration Code*	[] []
Concentration	[0] [28]
Vol. Std. Neat	[9.0] [9.0]
Vol. Std. Diluted	[0] [0]
Vol. Diluent	[0] [0]

*User Defined

Ensure that the assay parameters exactly match those listed above.

Axis-Shield Diagnostics Ltd.
 Luna Place
 The Technology Park, Dundee
 DD2 1XA, Scotland, UK
 Tel: +44 (0) 1382 422000
 Fax: +44 (0) 1382 422088
 Email: homocysteine@axis-shield.com

www.homocysteine.org.uk
 www.axis-shield.com

