(Distributer BECKMAN COULTER, samo za profesionalnu primjenu na AU platformama društva BECKMAN COULTER (AU400, AU480, AU680, AU5800, **DxC 500** i DxC 700 AU))

Axis-Shield Diagnostics Ltd.
The Technology Park
Dundee DD2 1XA
Ujedinjena Kraljevina
Tel: +44 (0) 1382 422000
Faks: +44 (0) 1382 422088

HRVATSKI:

NAMJENA

Liquid Stable (LS) 2-Part Homocysteine Reagent namijenjen je za *in vitro* kvantitativno određivanje ukupnog homocisteina u humanom serumu i plazmi. Proizvod može pomoći u postavljanju dijagnoze i liječenju pacijenata za koje se sumnja da imaju hiperhomocisteinemiju i homocistinuriju. Samo za profesionalnu uporabu.

UPOZORENJE: Uzorci pacijenta koji su na terapiji lijekovima koji uključuju S-adenozil metionin mogu pokazati netočno povišene razine homocisteina. Pacijenti koji uzimaju metotreksat, karbamazepin, fenitoin, dušikov oksid, antikoagulative ili 6-azauridin triacetat mogu imati povišene razine homocisteina zbog djelovanja tih tvari na metaboličke putove. Pogledajte poglavlje OGRANIČENJA ZA UPORABU u uputama u pakiranju priloženima uz ovaj test.

SAŽETAK I OBJAŠNJENJE TESTA

Homocistein (HCY) je aminokiselina koja sadrži tiol, a nastaje unutarstaničnom demetilacijom metionina. Homocistein se transportira u plazmu u kojoj, pretežno u oksidiranom stanju, cirkulira vezan za proteine plazme u obliku miješanog disulfida proteinskog homocisteina s albuminom (protein-SS-HCY).¹⁻⁵ Prisutne su manje količine reduciranog homocisteina i disulfidnog homocisteina (HCY-SS-HCY). Ukupni homocistein (tHCY) predstavlja zbroj svih vrsta homocisteina pronađenih u serumu ili plazmi (slobodnih i vezanih za proteine). Homocistein se metabolizira u cistein ili metionin. Pri transsulfuraciji vitamina B6 homocistein se nepovratno razlaže u cistein. Najveći dio homocisteina remetilacijom se pretvara u metionin, uglavnom s pomoću folata i enzima metionin-sintaze ovisnog o kobalaminu. Homocistein se nakuplja i izlučuje u krv kada dođe do ometanja ovih reakcija.^{3,5} Vrlo povišene koncentracije ukupnog homocisteina pronađene su u ispitanika s homocistinurijom, rijetkim genetskim poremećajem enzima koji su uključeni u metabolizam homocisteina. U pacijenata s homocistinurijom javljaju se mentalna retardacija, rana arterioskleroza te tromboembolije arterija i vena.^{2,6} Također su uočeni i drugi manje ozbiljni genetski defekti koji dovode do umjereno povišenih razina ukupnog homocisteina.⁷⁻⁹

Odnos između povišenih razina homocisteina i kardiovaskularne bolesti bio je predmet ispitivanja epidemioloških istraživanja. Meta-analizom 27 tih ispitivanja, uključujući više od 4000 pacijenata, procijenjeno je da je povećanje ukupnog homocisteina od 5 $\mu\text{mol/l}$ povezano s omjerom izgleda za koronarnu arterijsku bolest (CAD) od 1,6 (95 %-tni interval pouzdanosti [CI], od 1,4 do 1,7) za muškarce i 1,8 (95 %-tni CI od 1,3 do 1,9) za žene, a omjer izgleda za cerebrovaskularnu bolest bio je 1,5 (95 %-tni CI od 1,3 do 1,9). Rizik povezan s povećanjem od 5 $\mu\text{mol/l}$ ukupnog homocisteina bio je istovjetan riziku od povećanja kolesterola za 0,5 mmol/l (20 mg/dl). Također, ustanovljena je visoka korelacija s bolestima perifernih arterija.¹⁰

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* Vrijednosti podešene za rad u µmol

REFERENCIJE

1. McCully KS. Vascular Pathology of Homocysteinemia: Implications for the Pathogenesis of Arteriosclerosis. *Am J Pathol* 1969;56:111-122
2. Malinow MR. Plasma Homocyst(e)ine and Arterial Occlusive Diseases: A Mini-Review. *Clin Chem* 1995;41:173-176
3. Ueland PM. Homocysteine Species as Components of Plasma Redox Thiol Status. *Clin Chem* 1995;41:340-342
4. Perry IJ, Refsum H, Morris RW, et al. Prospective Study of Serum Total Homocysteine Concentration and Risk of Stroke in Middle-aged British Men. *The Lancet* 1995;346:1395-1398
5. Finkelstein JD. Methionine Metabolism in Mammals. *J Nutr Biochem* 1990;1:228-237
6. Mudd SH, Levy HL, Skovby F. Disorders of Transsulfuration. U: Scriver CR, Beaudet AL, Sly WS, et al., eds *The Metabolic and Molecular Basis of Inherited Disease*. New York: McGraw-Hill, 1995;1279-1327
7. Clarke R, Daly L, Robinson K, et al. Hyperhomocysteinemia: An Independent Risk Factor for Vascular Disease. *N Engl J Med* 1991;324:1149-1155
8. Deloughery TG, Evans A, Sadeghi A, et al. Common Mutation in Methylenetetrahydrofolate Reductase. *Circulation* 1996;94:3074-3078
9. Schmitz C, Lindpaintner K, Verhoef P, et al. Genetic Polymorphism of Methylenetetrahydrofolate Reductase and Myocardial Infarction. *Circulation* 1996;94:1812-1814
10. Boushey CJ, Beresford SAA, Omenn GS, et al. A Quantitative Assessment of Plasma Homocysteine as a Risk Factor for Vascular Disease. *JAMA* 1995;274:1049-1057
11. Cattaneo M. Hyperhomocysteinemia, Artherosclerosis and Thrombosis. *Thromb Haemost* 1999;81:165-176
12. Ridker PM, Manson JE, Buring JE, et al. Homocysteine and Risk of Cardiovascular Disease Among Postmenopausal Women. *JAMA* 1999;281:1817-1821
13. Bostom AG, Silbershatz H, Rosenberg IH, et al. Nonfasting Plasma Total Homocysteine Levels and All-Cause and Cardiovascular Disease Mortality in Elderly Framingham Men and Women. *Arch Intern Med* 1999;159:1077-1080
14. Guttormsen AB, Svarstad E, Ueland PM, et al. Elimination of Homocysteine from Plasma in Subjects with Endstage Renal Failure. *Irish J Med Sci* 1995;164 (Suppl. 15):8-9
15. Bostom AG, Lathrop L. Hyperhomocysteinemia in End-stage Renal Disease: Prevalence, Etiology, and Potential Relationship to Arteriosclerotic Outcomes. *Kidney Int* 1997;52:10-20
16. Ueland PM, Refsum H, Stabler SP, et al. Total Homocysteine in Plasma or Serum: Methods and Clinical Applications. *Clin Chem* 1993;39:1764-1779
17. Ueland PM, Refsum H. Plasma Homocysteine, A Risk Factor for Vascular Disease: Plasma Levels in Health, Disease, and Drug Therapy. *J Lab Clin Med* 1989;114:473-501
18. Fiskerstrand T, Refsum H, Kvalheim G, et al. Homocysteine and Other Thiols in Plasma and Urine: Automated Determination and Sample Stability. *Clin Chem* 1993;39:263-271
19. Nehler MR, Taylor LM Jr, Porter JM. Homocysteinemia as a Risk Factor for Atherosclerosis: A Review. *Cardiovascular Pathol* 1997;6:1-9
20. Lussier-Cacan S, Xhignesse M, Piolot A, et al. Plasma Total Homocysteine in Healthy Subjects: Sex-Specific Relation with Biological Traits. *Am J Clin Nutr* 1996;64:587-593
21. Clarke R, Woodhouse P, Ulvik A, et al. Variability and Determinants of Total Homocysteine Concentrations in Plasma in an Elderly Population. *Clin Chem* 1998;44:102-107
22. Jacques PF, Selhub J, Bostom AG, et al. The Effect of Folic Acid Fortification on Plasma Folate and Total Homocysteine Concentrations. *N Engl J Med* 1999;340:1449-1454
23. Lawrence JM, Pettiti DB, Watkins M and Umekubo MA. Trends in Serum Folate after Food Fortification. *The Lancet* 1999;354:915-916
24. Herrmann W, Schorr H, Obeid R, et al. Disturbed Homocysteine and Methionine Cycle Intermediates S-adenosylhomocysteine and S-adenosylmethionine are Related to Degree of Renal Insufficiency in Type 2 Diabetes. *Clin Chem* 2005;51:1-7
25. Obeid R, Kuhlmann MK, Kohler H, et al. Response of Homocysteine, Cystathionine, and Methylmalonic Acid to Vitamin Treatment in Dialysis Patients. *Clin Chem* 2005;51:196-201
26. Refsum H, Smith AD, Ueland PM, et al. Facts and recommendations about total homocysteine determinations: an expert opinion. *Clin Chem* 2004;50(1):3-32
27. National Committee for Clinical Laboratory Standards. *Method Comparison and Bias Estimation using Patient Samples; Approved Guideline-Second Edition*. NCCLS document EP9-A2. Wayne, PA: NCCLS, 2002
28. National Committee for Clinical Laboratory Standards. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2, Wayne, PA: NCCLS, 2004
29. National Committee for Clinical Laboratory Standards. *Protocols for the Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. NCCLS Document EP17-A. Wayne, PA: NCCLS, 2004.
30. Clinical Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP7-A2. Wayne, PA: CLSI, 2005.
31. Clinical Laboratory Standards Institute. *Measurement Procedure Comparison And Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP9-A3. Wayne, PA: CLSI, 2013.
32. Clinical Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Methods; Approved Guideline—Third Edition*. CLSI Document EP5- A3, Wayne, PA: CLSI, 2014
33. Clinical Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline - 2nd Edition*. CLSI Document EP17-A2, Wayne, PA: CLSI, 2012

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BC Distribution B.V.
Abbott Rapid Dx International Limited